STRUMENTARIO GALATA'

PER CHIRURGIA OSSEA

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The electrically driven modern instruments for bone diaeresis allow operations to be performed on the skeleton much quicker, with much more exactness and with much less traumatic effects than what could be done with the old instruments in the past.

For the purpose, of avoiding troubles that might have to weigh against the great advantages they offer, their construction must respect particular conditions however, to wit:

- 1) all the instruments and implements must have the surgical characteristics of lightness, size and shape, without losing the required strength;
- 2) their handling must be surgical and not mechanical, and their interchangeability must be quick;
- 3) the working speed must be adjustable by degrees, and the power — which must be widely sufficient for the resistances to be got over — must be constant under load;
- 4) the implements must be such as will allow the operator to stop and start their working, while the engine is operating, and he must have the direct drive of them;
- 5) the working must be entirely free from vibrations (which are so troublesome for operators) no matter what the speed may be;
- 6) sterilisation must be possible with steam, notwithstanding the presence of mechanical devices inside of the instruments;
- 7) and finally there must be absolute guarantees of safety, in spite of the usually high speed of the action.

The particular characteristics that make a distinction for this set of instruments are summed up at pag. 23, after en exhibit of the different parts.

A separate double-acting rotary pump, for suction of pathologic fluids and for lavages of cavities, can also be fitted to the engine of the instrument set by means of a back sled connecting piece.

LIST OF THE PARTS MAKING UP THE INSTRUMENT SET

Adjustable speed electric engine
Pedestal for the engine
Moulded metal case, containing the set of instruments
Flexible transmission shaft
Ordinary holder for the different implements
Angled holder for the circular saws
Automatic stop device trephine, with 2 borers and crown saw
Ossivorous screw forceps.

30 assorted implements and borers, made up as follows :

8 helicoidal pins for bone sutures and tacking (mm. 2, 2.5, 3, 3.5);
2 long helicoidal graduated borers (mm. 7 and 8);
6 spherical borers of Doyen (mm. 4, 6, 8, 10, 12, 14);
1 flat borer (mm. 12);
2 borers of Sudeck that can be graduated;
1 thin tubular saw;
1 trocar - trephine;
2 axles for circular saws;
1 double axle, that can be graduated, for coupled circular saws;
1 pair of surgeon's rasps (concave and

convex) for articular reconstructions;
2 other revolving rasps for moulding of the bones and cartilage;
2 revolving files.

22 circular saws, as follows :

2 sets of 8 saws (mm. 18, 20, 22, 24, 26, 30, 35, 40: hole of 6 mm.);
5 saws for coupled working (mm. 24, 26, 30, 35, 40: hole of 9 mm.);
1 saw of 45 mm. for cross cuts of large diaphyses.

The following accessories :

1 bent ribbon holder guide, with 2 protective metallic ribbons;
1 cut demarcation device;
1 elastic claw retractor;
2 bone levers of Tayel;
1 measurer of the thickness of the cranium;
1 wire holder;
2 supports to change the distance between the coupled saws, in the angled holder;
1 light mouse teeth pincers;
2 double spanners, 1 screw-driver, 1 flexometer, 1 oiler, 1 case with grease fusible at 150°, 2 metal brushes (hand and rotating) for bone splinters.

LIST OF THE PUMP'S OBJECTS

Rotary, double-acting pump, for suction and lavage
Sled connecting piece, to fit the pump to the engine
Pair of bottles, of 2 liters each, graduated, with large neck and large metallic closing

Metalic cage for the bottles, and accessories to apply it to the pedestal of the engine
Four metallic tubes (two angled), with detachable heads
Case containing the pump and the accessories.

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FIG. I ENGINE

Speed adjuster without rheostat (from 100 to 4000 revolutions per minute)
 Rod-insulator for power under load
 Automatic current switch
 Double repulsion junction (for the instrument set and for the pump)
 Closed frame

FIG. 2

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
FIG. 2

PEDESTAL

Height that can be regulated
Turning platform
Telescope support arm, for the
flexible transmission shaft


FIG. 3

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
FIG. 3

**ASSEMBLED ENGINE
AND FLEXIBLE
SHAFT**

(The support arm is to be used, during working, to hold the back part of the flexible shaft — when necessary and for the required length — and to hold the instrument after working).


FIG. 4 SCATOLA METALLICA
contenente lo strumentario

FIG. 4 MOULDED METAL CASE
containing the set of instruments


FIG. 5

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FIG. 5

COVER OF THE BOX
with implement-holder
frame, that can be detached

FIG. 6

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
FIG. 6

**AUTOMATIC STOP DEVICE
TREPHINE**

Automatic device (following figure).

Implements :

- 2 flat borers for craniotomy of 12 mm.;
- 1 crown saw of 18 mm. for endocranial decompressions and small explorations;
- 1 helicoidal point of 3 mm. for support sutures of cranial flaps (and to centre the crown).


FIG. 7

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
FIG. 7

Scheme of operation of the automatic trephine

During operation the stop ring goes up the instrument exactly to the same extent as the borer penetrates the bone. The borer stays motionless and the stop ring blocks further penetration of the borer at the same instant in which penetration of a bone of unknown thickness is finished. The automatic stop mechanism is instantaneous and it is absolutely sure mechanically. The depth of the maximum penetration of the borer is 20 mm.

The borer of 12 mm. can be quickly replaced by the crown saw of 18 mm., whose operation is also automatic.


FIG. 8


FIG. 9

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Application of one saw (fig. 8)

Application of coupled saws (fig. 9)

Converter set for the direction and level of the motion

Direct drive stop and starting device

Supports to change the distance between the coupled
circular saws

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
FIG. 11

How to fix and remove the saws on the angled holder

- 1) Keep the saw holder board motionless by keeping the front button of the instrument pressed downwards;
- 2) when the saw is fixed, screw up the nut by hand and, lastly, tighten moderately with the spanner.

To remove the saw, carry out the last named operation by starting at the end, but always keep the front button pressed downwards.

How to use the angled holder

- 1) Keep the fist low;
- 2) cut depthwise first;
- 3) only cut lengthwise afterwards and draw the instrument towards you.

FIG. 12

PUNTE ELICOIDALI

per suture e chiodature
ossee

(2 da 2 mm., 2 da 2,5 mm., 2
da 3 mm., 2 da 3,5 mm.)


FIG. 12

HELICOIDAL PINS

for bone sutures and
tacking

(2 of 2 mm., 2 of 2,5 mm., 2 of
3 mm., 2 of 3,5 mm.)

FIG. 13

FRESE SFERICHE E PIANA

(fresa sferiche da 4, 6, 8, 10, 12,
14 mm.; fresa piana da 12 mm.)


FIG. 13

SPHERICAL BORERS AND
FLAT BORER

(spherical borers of 4, 6, 8, 10,
12, 14 mm.; flat borer of. 12 mm.)

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FIG. 14 ORDINARY HOLDER
for all the implements

Direct drive stopping and starting device

FIG. 15



FIG. 16

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How to fit the implements on the ordinary holder (fig. 15). 1) Keep the instrument's lower button pressed backwards; 2) loosen the large front nut; 3) then run the borer into the holder up to the round groove cut on the stem and make it revolve slowly; 4) then tighten the nut by hand.

How to remove the implements from the ordinary holder (fig. 16). Keep the instrument's lower button pressed backwards 2) loosen the nut by hand and then follow out the foregoing operation backwards. Loosen with the spanner if the borer is tightly fixed after working.


FIG. 17

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AXLER FOR CIRCULAR SAWS
(to be mounted on the ordinary holder)

High resistance tungsten steel saws with special teeth to prevent pulpy deposits (saws of 18, 20, 22, 24, 26, 30, 35, 40, 45 mm.)


FIG. 18

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DOUBLE AXLE FOR COUPLED CIRCULAR SAWS

(to be mounted on the ordinary holder)

Sled pitch piece to be fixed by hand.
(Saws as above : of 24, 26, 30, 35, 40 mm.)


FIG. 19

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**BORERS SUDECK'S THAT CAN
BE GRADUATED**

for erosion cuts of flat bones
(2 of 4 mm.)

Jacket to regulate the depth (fixed, at
the required height, by the large nut
of the holder).


FIG. 20

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**PAIR OF SURGEON'S RASPS
(CONCAVE AND CONVEX)**

for articular reconstructions
(diarthrosis)


FIG. 21

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SURGEON'S RASPS AND FILES

for other moulding of the bones and
cartilage


FIG. 22

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**LONG AND GRADUATED
HELICOIDAL BORERS**

for infibulation of the neck of
the femur
(1 of 7 mm., 1 of 8 mm.)


FIG. 23 TREQUARTI - TRAPANO

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TROCAR - TREPINE

to take samples of bone medulla
(sternum, tibia and the large trochanter)


FIG. 24

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TUBULAR SAW

for taking samples of bones in
deep-seated places
Centering and expelling style


FIG. 25

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FIG. 25

OSSIVOROUS SCREW FORCEPS
for pressure cuts of flat bones


FIG. 26

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**RIBBON-HOLDER GUIDE AND PROTECTIVE
METALLIC RIBBON**

to protect the meninges during cuts of the sinciput

The front end of the ribbon-holder guide is slightly turned upwards and when it is in the entrance foramen of the sinciput, it makes the blunted bill of the flexible metallic ribbon likewise go upwards and get into immediate contact with the cranial vault up to the exit foramen.


FIG. 27

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CUT DEMARCACTION DEVICE

to demarcate the underlying protected sector exactly during cutting

When the ribbon is put in and the guide detached from it, the two ends (that come out of the foramina) are sent into a long metal handle, whose width is exactly the same as that of the underlying ribbon, so as to be able to demarcate the underlying protected sector exactly during the cutting.


FIG. 28

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Elastic claw retractors - Device to measure the thickness of
the cranium - Two Tavel levers - Wire holders.


FIG. 29

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Greasing under pressure for the instruments
(grease fusible at 150°)


FIG. 30

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RÉSUMÉ

OF THE CHARACTERISTICS OF THIS SET OF INSTRUMENTS

Originality of models

The engine, the four instruments and some implements and accessories are of original patented models.

Surgical measurements and shape of the instruments

The measurements, the weight and shape of the instruments have surgical characteristics.

Constant power - Adjustable speed

The power of the engine is just about constant under load, at any speed and with any implement. The engine runs with adjustable speed without a rheostat. The current is switched off automatically at zero.

Direct drive operation

The borers and saws can be stopped and started when the operator wishes, by raising or lowering a button with the thumb of the hand holding the instrument.

Safety stop device

In cases of need, the flexible transmission shaft can be expelled immediately from the working engine (and therefore reduced to immobility at once with the instrument and the implement) by merely raising a button.

No vibrations

The operation of any instrument or implement is free from vibrations no matter what the speed may be, on account of the exact centering of the various parts that go from the engine to the implement.

Speed and ease with cutting

Even on thick and compact bones the cuts of the circular saws are quick and easy, owing to the adequate power available, and the special teeth of the saws avoid pappy and splinter deposits forming between the saw and the bone, which lessen resistance considerably.

On account of this last named consideration, the saws do not get hot even after working for a long time at a high speed, and therefore necrosis due to burning is avoided.

The average cutting speed is 5 seconds per lineal centimeter of flat bone one centimeter thick. The cutting width is 1 mm.

Cutting in a perpendicular and axial direction

Cutting can take place not only in a direction perpendicular to the axis of the instrument, like what is done usually by means of the ordinary holder and the saw-holder axles, but also in an axial direction and at a lower level, with another special angled holder.

It is therefore possible to cut with the saw as one would cut with a bistoury and also execute, with perfect verticalness, transversal sections of the diaphysis and deep cuts.

Automatic stop perforations

A special trephine allows automatic stop perforations of the sinciput to be performed as well as of any flat bone of unknown thickness (not more than 20 mm.). The borer (or the crown) is blocked with strict automatism in the very same moment in which the perforation is finished.

N. B. - The aforesaid trephine, with a well studied ensemble of implements and accessories, allows craniotomy to be performed quickly, safely and with total elimination of the moments in which a trauma could take place.

Large equipment of instruments, implements and accessories

Besides the ordinary operations of bone separation, this equipment will allow particular or secondary bone separation (attic's cut, infibulation of the neck of the femur, endocranial decompressions and explorations, articular reconstructions of diarthrosis, moulding of bones and cartilage in general, to take deep-seated bone samples quickly and samples of bone marrow, etc.).

Steam sterilisation

The characteristics of the material allow sterilisation of the instruments even with flowing steam in spite of the mechanical devices inside of them.

Guarantee for lasting

The selection of the metals the strength of the construction and *suitable devices for inside lubrication with pressure* are a guarantee that Galata's set of instruments will last a long time and work well.

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FIG. 32

**PUMP FITTED ON THE ENGINE,
WITH THE ACCESSORIES**

Suction and washing can take turns, when required, through the two separate vessels by merely moving the drive lever sidewise.


FIG. 33

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FIG. 33

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FIG. 34 How to uncouple the pump from the engine

Raise the tenon button of the engine and the pump will be pushed backwards running on cylindrical guides.
It is sufficient to push the pump towards the engine to be able to couple.
In this way it is possible to obtain separately the operation of the pump or of the set of instruments, or operation, at the same time of the set of instruments and the pump together.


FIG. 35

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CASE CONTAINING THE PUMP AND THE ACCESSORIES

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NOTICE

1. An ordinary light socket will be sufficient to furnish electrical power to the Engine.
2. Detailed instructions for use will be sent with the set of Instruments and Pump.
3. Any special implement will be constructed on request, by the House.

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SOME OPERATIONS

FIG. 36
CROSS DIAPHYSIS SECTION

(For the great diaphysis:
cut on two sides)

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FIG. 37

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
FIG. 37

First stage

Parallel and longitudinal
cuts with the coupled saws
fixed on the ordinary
holder.

FIG. 38

Second stage

Cross cuts, at the two ends,
with the single saw fixed
on the angled holder.

FIG. 39

Third Stage

Removal with Tavel le-
vers of the isolated patch.

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
Fig. 40

First stage

Automatic perforations of the
sinciput

FIG. 41

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
FIG. 41

Second stage

- 1) Introduction of the protecting ribbon by means of the guide.

FIG. 42

- 2) Application of demarcation device on the protecting ribbon (after removing the guide).

FIG. 43

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
FIG. 43

3) Cuts of the bone patch on three sides.

FIG. 44

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
FIG. 44

Third stage

Cut of the end parts of the fourth side, by means of osteomyolic screw forceps, after retraction of the soft parts with the elastic claw retractor.

FIG. 45

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
FIG. 45

Fourth stage

- 1) Fracture of the middle part of the joint, with Tavel levers.

FIG. 46

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
FIG. 46

- 2) Turning over of the osteoplastic flap.

31455

55421